





























Ferritin ELISA

**Enzyme Immunoassay for Quantitative Determination of
Ferritin in Human Serum**

Instructions for use

1. SYMBOL LEGEND

	In vitro diagnostic medical device		EC Declaration of conformity
	Catalogue number		Batch code
	Use by		Manufacturer
	Date of manufacture		Consult operating instructions
	Temperature limitation		Biological risks
	Authorized Representative in the European Community		Conjugate
	Coated microplate (96 wells)		Substrate
	Wash solution, 20X concentrated		Stop solution
			
	Calibrators		Control
	Trial, 5000X concentrated		Optical density
			
	Contains sufficient for <n> tests		Deionized or distilled water
	Reconstitute with specified volume of liquid		Irritant
		Warning	

2. INTENDED USE

Ferritin kit is provided for the **quantitative determination of ferritin in human serum**.

Ferritin is an intracellular protein. It makes soluble iron pool that is necessary for successful erythropoiesis. At the same time ferritin protects the organism from the toxic effect of iron (connected with catalysis of free radicals production).

Ferritin molecule consists of iron-containing core and protein shell (apoferritin) with a molecular mass of about 450 kDa. Apoferritin consists of 24 subunits of “light” (L) and “heavy” (H) types. Each ferritin molecule can bind up to 4500 iron atoms in hydroxide and phosphate complexes.

Ferritin is located mainly in spleen, liver and red bone marrow blood corpuscles; in small quantities it is also present in plasma. Ferritin concentration in plasma usually permits to estimate total iron reserve in the organism adequately.

Ferritin concentration at the moment of birth is high (up to 600 ng/mL). In the first months of life it falls and remains at the level of about 30 ng/mL till pubescence. Then ferritin concentration slowly begins to grow, the level reached at the age of 24–25 remains for life. Normal concentration of ferritin in serum is 50–250 ng/mL for men and 20–150 ng/mL for women.

Ferritin concentrations below 10 ng/mL is a sign of asiderotic anemia. This parameter permits successfully differentiate asiderotic anemia from another types of anemia. Measurement of serum ferritin is also useful for monitoring of iron reserve in pregnant women, donors and patients subjected to regular haemodialysis. In the case of iron overloading ferritin concentration exceeds 400–500 ng/mL. It may reach several thousands ng/mL in the case of distinct hemochromathosis.

Concentration of serum ferritin increases in the case of infection, inflammatory processes (such as osteomyelitis and rheumatoid arthritis), acute and chronic liver disorders, leukemia, Hodgkin disease, breast cancer and some other oncopathologies.

3. PRINCIPLE OF THE TEST

Ferritin kit is a “sandwich” type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of ferritin molecule. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other coated onto on the inner surface of microwells. Ferritin molecules from the serum sample is bound to both immobilized antibody and anti-ferritin-peroxidase conjugate. Then the wells are washed with wash solution to remove any material not bound to the inner surface of the wells. Quantity of the bound to conjugate is directly proportional to ferritin level in the sample (Fig.1).

During the incubation with TMB solution the colour is developing. The intensity of the colour is directly proportional to ferritin concentration in specimens. Ferritin concentration in the patient sample is read from a standard curve that is processed in each assay




Fig. 1. Assay scheme

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label.

Ferritin kit should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 18 months.

After initial opening the kit is stable until the expiration date if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C for 12 months; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until the expiration date, protected from light;
- vials with concentrated Trial, concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8 °C for no more 4 weeks, in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle.

Damaged Test Kits

In case of any severe damage of the test kit or components, it has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 2 days. Aliquot and freeze samples for longer storage (–20 °C and lower). Avoid repeated freezing.

6. EXPECTED VALUES

Serum samples taken from 184 healthy men (age 21–45); 140 healthy women (age 19–43) and 52 pregnant women (age 20–34) were assayed with ferritin assay. The results are shown in the table below. These limits should be considered as guidelines only.

Study Group	Ferritin concentration range (ng/mL)	Mean ferritin concentration (ng/mL)
Healthy men	22-346	158
Healthy women	10-147	73
Pregnant women (I trimester)	55-90	61
Pregnant women (II trimester)	25-74	32
Pregnant women (III trimester)	10-16	12

It is highly recommended for each laboratory to determine its own reference range of ferritin concentrations.

7. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

8. REAGENT PREPARATION

Allow all the reagents to reach room temperature, and then thoroughly stir.

MP Keep **microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

CAL CONTROL Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature (+18...+25 °C) without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

WASH P Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL of **WASH P 20X** + 95 mL of water.

Mix thoroughly, avoid foaming.

SUB Protect **substrate** from direct light.

9. SAMPLE PREPARATION

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

10. PERFORMANCE CHARACTERISTICS OF THE ASSAY

10.1. Calibration - Traceability

Ferritin kit was calibrated against the WHO 2nd International Standard 80/578.

10.2. Specificity

According to Product specification provided by the supplier, no cross-interaction of both anti-ferritin monoclonal antibodies with serum albumin, haemoglobin, transferrin, alpha-fetoprotein and iron chloride was detected.

10.3. Analytical sensitivity

Analytical sensitivity of **Ferritin, kit** i.e. concentration that can be distinguished from zero calibrator, is 5 ng/mL. It is defined as mean OD of 10 replicates of Calibrator 0 plus two standard deviations.

10.4. Measurement Range

Ferritin kit was validated for measurement of ferritin concentration within the concentration diapason of 5-1000 ng/mL.

10.5. Hook Effect

For **Ferritin kit high dose hook effect** was not detected for concentrations up to 10 000 ng/mL. **High dose hook effect** was determined by spiking calibrator 0 matrix with antigen.

10.6. Intra- and Inter-Assay Variation

To determine **intra-assay CV** 8 serum samples were assayed in 9 replicates each. The results are shown below.

Sample	Mean ferritin concentration, ng/mL	Intra-assay CV	
		SD	CV, %
HS 1	16.4	0.62	3.8
HS 2	18.1	0.64	3.5
HS 3	59.9	1.85	3.1
HS 4	113	2.4	2.1
HS 5	151	3.7	2.4
HS 6	343	21.3	6.2
HS 7	450	19.8	4.4
HS 8	627	12.5	2.0

To determine **inter-assay CV** 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are shown below.

Sample	Mean ferritin concentration, ng/mL			Inter-assay CV	
	1 assay	2 assay	3 assay	SD	CV, %
HS 1	15.2	13.7	14.1	0.79	5.5
HS 2	44.0	44.8	43.2	0.90	2.0
HS 3	51.3	49.1	54.2	2.52	5.0
HS 4	81.0	80.5	80.1	0.50	0.6
HS 5	120	124	120	2.1	1.8
HS 6	277	275	283	4.2	1.5
HS 7	366	375	382	8.0	2.2
HS 8	583	588	609	13.8	2.4

11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

12. SAFETY PRECAUTIONS

• **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.

• Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.

• Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.

• Do not use substrate, stop solution and wash solution supplied by other vendors.

• Use only “P”-labeled wash solution.

• Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.

• Take into account the following common procedural notes:

-always pipette reagents into wells immediately after washing procedure;

-avoid contamination of the solutions;

-in case of partial use of the kit, dispense only required volume of the reagent into the tray; do not pour unused reagents back into the original vials;


-avoid exposure to direct sunlight during incubations;


-always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;


-the incubation temperature for all the immunological reactions must be kept at +37 °C;


-do not touch the bottom of the wells;
-calibrators should be measured in each separate assay. It is also recommended to measure each time ferritin concentration in the control.

- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.

-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.

-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.

-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.

-  As the kit contains potentially hazardous material, the following precautions should be taken:

- do not smoke, eat or drink while performing the assay;

- always use protective gloves;
- never pipette material by mouth;
- in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.

• GLP including all general and individual regulations should be applied for the kit usage.



As the kit contains irritant (**CONJ**, **CAL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

13. MATERIAL PROVIDED

13.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-ferritin monoclonal antibodies.	1 pcs.
CONJ	Conjugate: solution contains anti-ferritin monoclonal antibodies conjugated with HRP.	14 mL, ready to use
0-5 CAL	Ferritin calibrators: protein-based solution or lyophilized preparations containing known ferritin concentrations – 0;10;30;100;300;1000 ng/mL (approximate values). The concentrations of calibrators may be different for schemes with or without shaking. For exact ferritin concentrations, see vial labels.	6 vials, 0.5 mL each ready to use or lyophilized preparations.
CONTROL	Ferritin control: protein-based solution or lyophilized preparation containing known ferritin concentration. The range of ferritin concentration may be different for schemes with or without shaking. For exact range of ferritin concentration, see vial label.	0.5 mL, ready to use or lyophilized preparation
SUB	Substrate (TMB solution): 3, 3', 5, 5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide.	14 mL, ready to use
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 560 of solution.	2x14 mL, concentrated
STOP	Stop solution: 1 N HCl solution.	14 mL, ready to use

13.2. Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator/shaker (+37 °C, shaking speed 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

13.3. Test Procedure

Ferritin kit is designed for 96 tests. This is sufficient for 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

13.3.1. Assay Procedure

13.3.1.1. Protocol with shaking

(See assay scheme, section 13.5.)

All samples should be tested in duplicates.

A. Pipette:

- 20 μL of calibrators **CAL**, control **CONTROL** and patient's samples in duplicates into the respective wells;

Leave wells A1-A2 empty for blank!

B. Pipette:

- 100 μL of conjugate **CONJ** into each well, except wells A1-A2

Note: total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.

C. Incubate for 30 minutes at +37 °C while shaking (500–800 rpm).

D. Wash 5 times, as described in section 13.3.2.

E. Pipette 100 μL of substrate **SUB into each well (including blank); incubate strips at room temperature (+18...+25 °C) in the dark for 15-30 minutes, depending on the colour intensity, or 10 minutes while shaking (500-800 rpm) at +37 °C.**

F. Pipette 100 μL of stop solution **STOP into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. Shake for 1–2 min at room temperature.**

G. Read OD at 450 nm within 20 min.

13.3.1.2. Protocol without shaking

(See assay scheme, section 13.6.)

All samples should be tested in duplicates.

A. Pipette:

- **20 µL** of calibrators **CAL**, control **CONTROL** and patient's samples in duplicates into the respective wells;

Leave wells A1-A2 empty for blank!

B. Pipette:

- **100 µL** of conjugate **CONJ** into each well, except wells A1-A2

Note: total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.

C. Incubate for 60 minutes at 37°C (pre-shake 1-2 minutes at room temperature).

D. Wash 5 times, as described in section 13.3.2.

E. Pipette 100 µL of substrate SUB into each well (including blank); incubate strips at room temperature (18-25°C) **in the dark for 15-30 minutes**, depending on the colour intensity.

F. Pipette 100 µL of stop solution STOP into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. Shake for 1–2 min at room temperature.

G. Read OD at 450 nm within 20 min.

13.3.2. Wash Procedure

It is advisable to use an automatic microplate washer set at 5 wash cycles and a volume of 300 μL of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;
- dispense 300 μL of wash solution, prepared according to section 8, into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 5 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

13.4. Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated = $2.28 - 0.06 = 2.22$

13.4.1. Data Reliability (for OD Measured at 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) ≤ 0.100 ;
- average OD of Cal 5 ≥ 1.5 (after blank subtraction).
- control's concentration must fall within the acceptability range that is shown on the vial label.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

13.4.2. Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective ferritin concentrations using 4PL fit (see typical standard curve, fig. 2). Calculate concentration of ferritin in samples using standard curve.

Any extrapolation of the standard curve to ferritin above the nominal value of the calibrator 5 (approximately 1000 ng/mL) is forbidden.




Fig. 2. Example of typical standard curve.
Do not use for evaluation of real assay data!

13.5. Assay scheme with shaking

	Wells		CAL	
Reagents		«Blank»	CONTROL	Samples
CAL		–	20 µL	–
CONTROL		–	–	–
Samples		–	–	20 µL
CONJ		–	100 µL	100 µL
Incubation No.1	30 min, +37 °C, 500–800 rpm			
WASH P (diluted)	5 x 300 µL			
SUB	100 µL	100 µL	100 µL	100 µL
Incubation No.2	15–30 min, +18...+25 °C, in the dark			
	10 min, +37 °C, 650-800 rpm			
STOP	100 µL	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C			
OD measuring	450 nm			
Calculations	Corresponding software			

13.6. Assay scheme without shaking

Wells	«Blank»	CAL CONTROL	Samples
Reagents			
CAL CONTROL	–	20 µL	–
Samples	–	–	20 µL
CONJ	–	100 µL	100 µL
Incubation No.1	60 min, 37°C (pre-shake 1-2 minutes at room temperature)		
WASH P (diluted)	5 x 300 µL		
SUB	100 µL	100 µL	100 µL
Incubation No.2	15–30 min, 18-25°C, in the dark		
STOP	100 µL	100 µL	100 µL
Stirring	1–2 min, 18-25°C		
OD measuring	450 nm		
Calculations	Corresponding software		

BIOGENIX INC. PVT. LTD.

Factory: B - 19/A, S.I.L Ancillary Estate, Amausi Industrial Area, Nadarganj, Kanpur Road, Lucknow - 226008 (U.P.), India

Email: biogenix2007@yahoo.com, info@biogenixinc.com

Website: www.biogenixinc.com

Customer care no: +919140971443