



## SINGLE CHANNEL TYPE FLUORESCENT IMMUNO ANALYZER (POCT)

BI-180

### BASIC PARAMETERS

<b>Weight</b>	2.5 kg
<b>Size</b>	280 mm (L) * 220 mm (W) * 150 mm (H)
<b>Screen</b>	7.05" 24-bit true color touch-screen
<b>Speed</b>	180 T/h
<b>Sample Types</b>	serum/plasma, whole blood, urine, etc.
<b>Data Transmission</b>	USB, RS-232, wifi, ethernet, support LIS and HIS
<b>Print</b>	built-in thermal printer, can attach to external printer
<b>Storage</b>	> 10,000
<b>Operation</b>	automatically test and discard reagent cards; automatically identify project ID and information;
<b>Quality Control</b>	double quality control of instrument and reagents